



PRIZJet



ISO 9001 Q-1111
ISO 13485 M-017
ISO 14001 E-397



COMPRESSOR

NEBULISER

COMPRESSOR

NEBULISER

- QUALITY AEROSOL
- NEB KIT HOLDER
- INTERMITTENT OR CONTINUOUS MODE OF DRUG DELIVERY
- CARRY BAG
- ERGONOMIC DESIGN
- EASY TO USE

CE 1304



The **PRIZ Jet** is a high quality, mains operated compressor nebulizer. It is easy to use for reliable respiratory disorder treatment at home or in the hospital. Unique neb kit provides two modes of drug delivery – intermittent or continuous. The **PRIZ Jet** is equipped with a set of accessories needed for adult and children highly effective treatment.

PRIZ Jet with accessories



Specification:

- Product name: PRIZMA compressor nebuliser
- Model: PRIZ Jet™
- Compressor maximum pressure: 2,5 bar
- Compressor air flow: 8 l/min
- Aerosol output: 0.3 ml/min
- Nebuliser capacity: 6 ml
- Residual volume: 0.64 ml
- Particle size MMAD: 3 µm
- Sound level: < 60 dB (A)
- Operating voltage: 230V/50 Hz
- Power consumption: 45 W
- Dimensions: 180 x 140 x 110 mm
- Weight: 1,3 kg (compressor only)