



PROFI III[®]
SONIC



ISO 9001 Q-1111
ISO 13485 M-017
ISO 14001 E-397



PROFESSIONAL
ULTRASONIC
NEBULISER

**PROFESSIONAL
ULTRASONIC
NEBULISER**

PROFESSIONAL ULTRASONIC NEBULISER

CE 1304



PROFI SONIC is ultrasonic nebuliser which is perfect for use in hospitals, nursing homes and medical clinics.

This high volume nebuliser with 1-100 minutes timer and large touch screen display makes it unique of its kind.

It is supplied with bacterial/viral filter offering wide range of accessories such as bottle for continual nebulisation, heatable aerosol hose, etc.

Touch screen display



- PERFECT FOR USE IN HOSPITALS, NURSING HOMES AND MEDICAL CLINICS
- HIGH VOLUME NEBULISER
- TIMER: 1 – 100 MINUTES
- LARGE TOUCH SCREEN DISPLAY
- SUPPLIED WITH ASSEMBLED BACTERIAL/VIRAL FILTER
- WIDE RANGE OF ACCESSORIES AVAILABLE (bottle for continual nebulisation, heatable aerosol hose)

Specification:

- Product name: PRIZMA ultrasonic nebuliser
- Model: PROFI SONIC™
- Ultrasonic frequency: 1,7 MHz
- Nebulisation rate: 0,0 – 3,0 ml/min.
- Particle size range: 0,5 – 5,0 µm
- Air volume: max 20 l/min.
- Medication cup capacity: 150 ml
- Water tank capacity: 400 ml
- Sound level: 35 dB (A)
- Operating voltage: 100-240 V AC, 50/60 Hz
- Power consumption: 50 W
- Dimension: 250 x 200 x 250 mm
- Weight: 2 kg